

Document made available under the Patent Cooperation Treaty (PCT)

International application number: PCT/US05/005597

International filing date: 24 February 2005 (24.02.2005)

Document type: Certified copy of priority document

Document details: Country/Office: US
Number: 60/547,168
Filing date: 25 February 2004 (25.02.2004)

Date of receipt at the International Bureau: 31 March 2005 (31.03.2005)

Remark: Priority document submitted or transmitted to the International Bureau in compliance with Rule 17.1(a) or (b)



World Intellectual Property Organization (WIPO) - Geneva, Switzerland
Organisation Mondiale de la Propriété Intellectuelle (OMPI) - Genève, Suisse

1298854

THE UNITED STATES OF AMERICA

TO ALL TO WHOM THESE PRESENTS SHALL COME:

UNITED STATES DEPARTMENT OF COMMERCE

United States Patent and Trademark Office

March 21, 2005

THIS IS TO CERTIFY THAT ANNEXED HERETO IS A TRUE COPY FROM THE RECORDS OF THE UNITED STATES PATENT AND TRADEMARK OFFICE OF THOSE PAPERS OF THE BELOW IDENTIFIED PATENT APPLICATION THAT MET THE REQUIREMENTS TO BE GRANTED A FILING DATE.

APPLICATION NUMBER: 60/547,168

FILING DATE: February 25, 2004

RELATED PCT APPLICATION NUMBER: PCT/US05/05597



Certified by

Under Secretary of Commerce
for Intellectual Property
and Director of the United States
Patent and Trademark Office

14023 U.S. PTO
022504

PROVISIONAL APPLICATION FOR PATENT COVER SHEET

This is a request for filing a PROVISIONAL APPLICATION FOR PATENT under 37 CFR 1.53(c).

Express Mail Label No. ET 889956556 US

INVENTOR(S)

Given Name (first and middle [if any])	Family Name or Surname	Residence (City and either State or Foreign Country)
Darron K. Michael L. Hong X.	Young Roukes Tang	South Pasadena, CA Pasadena, CA Pasadena, CA

15535 U.S. PTO
60/547168

022504

☐ Additional inventors are being named on the _____ separately numbered sheets attached hereto

TITLE OF THE INVENTION (280 characters max)

Wafer Scale Fabrication of Nanoscale Scanning Probes for Biological and High Frequency Applications

Direct all correspondence to :

CORRESPONDENCE ADDRESS

☒ Customer Number

000029690

OR

Place this label in the
Bar Code box here

29690

ENCLOSED APPLICATION PARTS (check all that apply)

PATENT TRADEMARK OFFICE

☒ Specification Number of Pages 4

☐ CD(s), Number

☐ Drawing(s) Number of Sheets _____

☐ Application Data Sheet. See 37 CFR 1.76

☒ Other (specify)

Postcard

METHOD OF PAYMENT OF FILING FEES FOR THIS PROVISIONAL APPLICATION FOR PATENT

☒ Applicant claims small entity status. See 37 CFR 1.27.

☐ A check or money order is enclosed to cover the filing fees.

☒ The Commissioner is hereby authorized to charge filing fees or credit any overpayment to Deposit Account Number:

50-1742

FILING FEE
AMOUNT (\$)

80

☐ Payment by credit card. Form PTO-2038 is attached.

The invention was made by an agency of the United States Government or under a contract with an agency of the United States Government.

☐ No

☒ Yes, the name of the U.S. Government agency and the contract number are: **Grant No. awarded by the DARPA**

Respectfully submitted,

Date 2/24/2004

SIGNATURE

Rich Wolf

REGISTRATION NO.

47,319

TYPED or PRINTED NAME: **Rich Wolf**

(if appropriate)

Docket Number:

CIT-4049-P

TELEPHONE: (626) 395-2322

Wafer scale fabrication of nanoscale scanning probes for biological and high frequency applications

Abstract:


The present art describes the design, fabrication and characterization of nanoscale cantilevers that can be used for scanning probe microscopy of biological samples. The dimensions of the cantilevers fall into nanoscale regime. For example, the thickness of the cantilever can be as thin as 30nm. The width of the cantilevers is as narrow as 400nm. These dimensions match the size of cells and offer high resolution in scanning probe microscopes (SPM). Meanwhile, the small dimensions of the cantilevers allow them to be operated at much higher frequencies and yield faster temporal response (as fast as 1uS).

The fabrication process employs a wafer scale nanofabrication technique that is capable of producing hundreds of cantilevers from one wafer. The silicon chips that house the cantilever match the size of the chips used in the SPM market. They also have electronic circuits that can interface directly into commercial scanning probe microscopes and provide force-sensitive electrical signal.

Two types of force transduction schemes are proposed and tested on these nanoscale cantilevers: piezoresistive detection and optical detection. The former relies on the high piezoresistivity of epitaxially grown p+ silicon layer embedded in the cantilever. Optical detection is realized by extending the nanocantilever with a larger micron-scale reflecting pads. Both methods have demonstrated high sensitivity in the fabricated devices

Material for piezoresistive detection:

The structure is 30 nm p-doped ($4 \times 10^{19} \text{ Be/cm}^2$) Si epilayer grown on 80nm undoped Si layer, on top of 400nm of SiO_2 , on 300 μm thick handle wafer oriented in [100] direction. This will maximize the piezoresistive constant along the [110] direction and ease fabrication and alignment of scanning probes.¹ See Figure 1.



Material for optical detection:


A similar structure is designed for optimal detection. Here the p-doped layer is eliminated.

Probe design for piezoresistive detection:


Align probe pattern along direction with maximum piezoresistive coefficient. Design scanning probe with electrode “legs” to form a stress concentration region where the current passes the most compliant area of the probe, maximizing piezoresistive response with actuation.

Probe design for optical detection:


Design scanning probe with at least $\sim 1\mu\text{m}^2$ reflective area near the tip of probe. Evaporate/sputter metal (i.e. 20nm of Au) on area near probe tip. This allows diffraction limited optical feedback to be utilized in commercial scanning probe microscopes (see Figure 2).



Release Structure Design:



Define temporary support structure for probe chip during nanoscale fabrication and two release areas enabling tweezer removal of chip. (See Figure 3). Two support beams are placed near the tip of chip to protect membrane from tearing prior to defining scanning probe tip. One support beam is use to support rear end of chip to wafer. Support beams have to be strong enough to withstand future wafer-scale processing of scanning probe while compliant enough to be broken for chip removal from wafer (i.e. 25 X 70 X .1 microns).




Process Flow:

- 1) Electrode fabrication: ohmic metal evaporation, define electrode pattern with photolithography, wet etch
- 2) Membrane and release-structure fabrication: Deep Reactive Ion Etching through Si handle wafer
- 3) Oxide removal: buffered oxide wet etch
- 4) Cantilever definition: electron beam/ photo lithography of mask
- 5) Cantilever fabrication: etch cantilever shape in piezoresistive membrane
- 6) Release of scanning probe chip: pull out chip from wafer with tweezers

BEST AVAILABLE COPY

Characterization:

Optically detected signal (Figure 4)



Piezoresistive detection signal (Figure 5)

